



Zimmer® Periarticular Proximal Tibial Locking Plates



Plate Features:

- Anatomically contoured plates create a fit that requires little or no additional bending and are helpful for metaphyseal/diaphyseal reduction
- The anatomical shape of the head of the plate matches the shape of the proximal tibia
- Locking holes in the plate head allow placement of the screws to capture medial fragments
- The low profile plate facilitates fixation without impinging on soft tissues
- Thick-to-thin profiles make the plates autocontourable
- The tapered plate shaft design allows for a minimally invasive technique with submuscular passage of the plate
- The plate can be used to support a medial wedge or fragment
- 3.5mm Proximal Lateral Tibial Locking Plates are available in six lengths, from 6 holes (104mm) to 16 holes (224mm)
- 5.5mm Proximal Lateral Tibial Locking Plates are available in six lengths, from 4 holes (97mm) to 14 holes (250mm)



Prod. No. (Left)	Prod. No. (Right)	Description	Holes	Length
3.5mm Proximal Lateral Tibial Locking Plates				
00-2357-004-06	00-2357-003-06	3.5mm Proximal Lateral Tibial Locking Plate	6	104mm
00-2357-004-08	00-2357-003-08	3.5mm Proximal Lateral Tibial Locking Plate	8	128mm
00-2357-004-10	00-2357-003-10	3.5mm Proximal Lateral Tibial Locking Plate	10	152mm
00-2357-004-12	00-2357-003-12	3.5mm Proximal Lateral Tibial Locking Plate	12	176mm
00-2357-004-14	00-2357-003-14	3.5mm Proximal Lateral Tibial Locking Plate	14	200mm
00-2357-004-16	00-2357-003-16	3.5mm Proximal Lateral Tibial Locking Plate	16	224mm

Designed for use with:

- Zimmer 3.5mm/2.7mm Periarticular Locking Screw Set (00-2359-000-02)
- Zimmer 3.5mm/2.7mm Periarticular Locking Instrument Set (00-2360-000-02)
- Zimmer 3.5mm Periarticular Proximal Lateral Tibial Plate Standard Jig Set (00-2360-000-15)

Compatible with:

- Zimmer Periarticular 3.5mm Self-tapping Cortical Screws
- Zimmer Periarticular 4.0mm Cancellous Screws
- Zimmer Small Fragment Instrument Set (00-4800-035-20)
- Zimmer Small Fragment Screw Set (00-4800-035-01)



Prod. No. (Left)	Prod. No. (Right)	Description	Holes	Length
5.5mm Proximal Lateral Tibial Locking Plates				
00-2357-006-04	00-2357-005-04	5.5mm Proximal Lateral Tibial Locking Plate	4	97mm
00-2357-006-06	00-2357-005-06	5.5mm Proximal Lateral Tibial Locking Plate	6	128mm
00-2357-006-08	00-2357-005-08	5.5mm Proximal Lateral Tibial Locking Plate	8	158mm
00-2357-006-10	00-2357-005-10	5.5mm Proximal Lateral Tibial Locking Plate	10	189mm
00-2357-006-12	00-2357-005-12	5.5mm Proximal Lateral Tibial Locking Plate	12	219mm
00-2357-006-14	00-2357-005-14	5.5mm Proximal Lateral Tibial Locking Plate	14	250mm

Designed for use with:


- Zimmer 5.5mm/4.5mm Periarticular Locking Screw Set (00-2359-000-01)
- Zimmer 5.5mm/4.5mm Periarticular Locking Instrument Set (00-2360-000-01)
- Zimmer 5.5mm Periarticular Proximal Lateral Tibial Plate Standard Jig Set (00-2360-000-12)

Compatible with:

- Zimmer Periarticular 4.5mm Self-tapping Cortical Screws
- Zimmer Periarticular 6.5mm Cancellous Screws
- Zimmer Basic Instrument Set (00-4800-045-20)
- Zimmer Standard Screw Set (00-4800-045-01)

Prod. No.	Description	Set Qty.
00-2357-000-09	3.5mm Proximal Lateral Tibial Locking Plate Set	
00-2357-004-06	3.5mm Proximal Lateral Tibial Locking Plate, 6 Hole, Left	1
00-2357-004-08	3.5mm Proximal Lateral Tibial Locking Plate, 8 Hole, Left	1
00-2357-004-10	3.5mm Proximal Lateral Tibial Locking Plate, 10 Hole, Left	1
00-2357-004-12	3.5mm Proximal Lateral Tibial Locking Plate, 12 Hole, Left	1
00-2357-004-14	3.5mm Proximal Lateral Tibial Locking Plate, 14 Hole, Left	1
00-2357-004-16	3.5mm Proximal Lateral Tibial Locking Plate, 16 Hole, Left	1
00-2357-003-06	3.5mm Proximal Lateral Tibial Locking Plate, 6 Hole, Right	1
00-2357-003-08	3.5mm Proximal Lateral Tibial Locking Plate, 8 Hole, Right	1
00-2357-003-10	3.5mm Proximal Lateral Tibial Locking Plate, 10 Hole, Right	1
00-2357-003-12	3.5mm Proximal Lateral Tibial Locking Plate, 12 Hole, Right	1
00-2357-003-14	3.5mm Proximal Lateral Tibial Locking Plate, 14 Hole, Right	1
00-2357-003-16	3.5mm Proximal Lateral Tibial Locking Plate, 16 Hole, Right	1

Prod. No.	Description	Set Qty.
00-2357-000-08	5.5mm Proximal Lateral Tibial Locking Plate Set	
00-2357-006-04	5.5mm Proximal Lateral Tibial Locking Plate, 4 Hole, Left	1
00-2357-006-06	5.5mm Proximal Lateral Tibial Locking Plate, 6 Hole, Left	1
00-2357-006-08	5.5mm Proximal Lateral Tibial Locking Plate, 8 Hole, Left	1
00-2357-006-10	5.5mm Proximal Lateral Tibial Locking Plate, 10 Hole, Left	1
00-2357-006-12	5.5mm Proximal Lateral Tibial Locking Plate, 12 Hole, Left	1
00-2357-006-14	5.5mm Proximal Lateral Tibial Locking Plate, 14 Hole, Left	1
00-2357-005-04	5.5mm Proximal Lateral Tibial Locking Plate, 4 Hole, Right	1
00-2357-005-06	5.5mm Proximal Lateral Tibial Locking Plate, 6 Hole, Right	1
00-2357-005-08	5.5mm Proximal Lateral Tibial Locking Plate, 8 Hole, Right	1
00-2357-005-10	5.5mm Proximal Lateral Tibial Locking Plate, 10 Hole, Right	1
00-2357-005-12	5.5mm Proximal Lateral Tibial Locking Plate, 12 Hole, Right	1
00-2357-005-14	5.5mm Proximal Lateral Tibial Locking Plate, 14 Hole, Right	1



Prod. No.	Description	Set Qty.
00-2360-000-15	3.5mm Proximal Tibial Standard Jig Set	
00-2360-093-01	3.5mm Proximal Lateral Tibia Plate Jig, Right	1
00-2360-093-02	3.5mm Proximal Lateral Tibia Plate Jig, Left	1
00-2358-050-00	Case	1

Prod. No.	Description	Set Qty.
00-2360-000-12	5.5mm Proximal Tibial Standard Jig Set	
00-2360-091-01	5.5mm Proximal Lateral Tibia Plate Jig, Right	1
00-2360-091-02	5.5mm Proximal Lateral Tibia Plate Jig, Left	1
00-2358-015-00	Case	1

Prod. No.	Description	Set Qty.
00-2360-000-02	3.5mm/2.7mm Periarticular Locking Instrument Set	
00-1147-071-00	Cleaning Stylet	1
00-1147-076-00	Cleaning Brush	1
00-2358-040-00	3.5/2.7mm Locking Screw and Instrument Case	1
00-2360-011-01	3.5mm/2.7mm Locking Plate Reduction Instrument	2
00-2360-011-02	3.5mm/2.7mm Locking Plate Reduction Sleeve	2
00-2360-012-03	Plate Reduction Spin Knob	2
00-2360-020-20	2.0mm Standard Cannula	4
00-2360-020-27	2.7mm Standard Cannula	4
00-2360-021-16	1.6mm Standard Cannula	4
00-2360-033-16	1.6mm Standard Drill Tip Guide Wire (5 per box)	2 bx
00-2360-040-35	3.5mm/2.7mm Locking Screw Standard Depth Gauge	1
00-2360-041-35	3.5mm Cannulated Locking Screw Depth Gauge	1
00-2360-053-27	2.7mm Locking Screw Tap	2
00-2360-053-35	3.5mm Locking Screw Tap	2
00-2360-054-35	3.5mm Cannulated Locking Screw Tap	2
00-2360-065-00	2.5mm Screwdriver Stop Ring	1
00-2360-065-25	2.5mm Hex Standard Screwdriver	2
00-2360-066-25	2.5mm Hex Standard Cannulated Screwdriver	2
00-2360-071-27	2.7mm Standard Cannulated Drill	2
00-2360-080-00	Torque Limiting Attachment	1
00-2360-086-00	Modular Handle	1
00-2360-088-00	3.5mm/2.7mm Cannula Insertor	1
00-2360-175-20	2.0mm Standard Drill	2
00-2360-205-27	2.7mm Standard Drill	2
00-4812-045-00	Large Hex Screwdriver	1
Also Available:		
00-2360-050-25	2.5mm Hex Forward Captive Screwdriver	0
00-2360-051-25	2.5mm Hex Cannulated Forward Captive Screwdriver	0

Prod. No.	Description	Set Qty.
00-2360-000-01	5.5mm/4.5mm Periarticular Locking Instrument Set	
00-1147-073-00	Cleaning Stylet	1
00-1147-078-00	Cleaning Brush	1
00-2358-035-00	5.5mm/4.5mm Locking Screw and Instrument Case	1
00-2360-012-01	5.5mm/4.5mm Plate Reduction Instrument	2
00-2360-012-02	5.5mm/4.5mm Plate Reduction Sleeve	2
00-2360-012-03	Plate Reduction Spin Knob	2
00-2360-020-37	3.7mm Standard Cannula	4
00-2360-021-32	3.2mm Standard Cannula	4
00-2360-033-32	3.2mm Standard Drill Tip Guide Wire (5 per box)	2 bx
00-2360-040-45	4.5mm Locking Screw Standard Depth Gauge	1
00-2360-041-55	5.5mm Cannulated Locking Screw Depth Gauge	1
00-2360-053-45	4.5mm Locking Screw Tap	2
00-2360-054-55	5.5mm Cannulated Locking Screw Tap	2
00-2360-065-05	5.0mm Screwdriver Stop Ring	1 bx
00-2360-065-50	5.0mm Hex Standard Screwdriver	2
00-2360-066-50	5.0mm Hex Standard Cannulated Screwdriver	2
00-2360-071-47	4.7mm Standard Cannulated Drill	2
00-2360-080-00	Torque Limiting Attachment	1
00-2360-085-00	Guide Wire Insertor	1
00-2360-086-00	Modular Handle	1
00-2360-087-00	Ratcheting Modular Handle	1
00-2360-088-00	5.5mm/4.5mm Cannula Insertor	1
00-2360-090-04	5.5mm/4.5mm Standard Jig Sleeve	4
00-2360-225-37	3.7mm Standard Drill	2
Also Available:		
00-2360-050-50	5.0mm Hex Forward Captive Screwdriver	0
00-2360-051-50	5.0mm Hex Cannulated Forward Captive Screwdriver	0

Prod. No.	Description	Set Qty.
00-2359-000-02	3.5mm/2.7mm Locking Screw Set	
00-2359-030-36	3.5mm Cannulated Locking Screw 30mm Lng	3
00-2359-035-36	3.5mm Cannulated Locking Screw 35mm Lng	3
00-2359-040-36	3.5mm Cannulated Locking Screw 40mm Lng	3
00-2359-045-36	3.5mm Cannulated Locking Screw 45mm Lng	3
00-2359-050-36	3.5mm Cannulated Locking Screw 50mm Lng	4
00-2359-055-36	3.5mm Cannulated Locking Screw 55mm Lng	4
00-2359-060-36	3.5mm Cannulated Locking Screw 60mm Lng	4
00-2359-065-36	3.5mm Cannulated Locking Screw 65mm Lng	4
00-2359-070-36	3.5mm Cannulated Locking Screw 70mm Lng	4
00-2359-075-36	3.5mm Cannulated Locking Screw 75mm Lng	3
00-2359-080-36	3.5mm Cannulated Locking Screw 80mm Lng	3
00-2359-085-36	3.5mm Cannulated Locking Screw 85mm Lng	3
00-2359-090-36	3.5mm Cannulated Locking Screw 90mm Lng	3
00-2359-030-37	3.5mm Cannulated Conical Screw 30mm Lng	2
00-2359-035-37	3.5mm Cannulated Conical Screw 35mm Lng	2
00-2359-040-37	3.5mm Cannulated Conical Screw 40mm Lng	2
00-2359-045-37	3.5mm Cannulated Conical Screw 45mm Lng	2
00-2359-050-37	3.5mm Cannulated Conical Screw 50mm Lng	2
00-2359-055-37	3.5mm Cannulated Conical Screw 55mm Lng	2
00-2359-060-37	3.5mm Cannulated Conical Screw 60mm Lng	2
00-2359-065-37	3.5mm Cannulated Conical Screw 65mm Lng	2
00-2359-070-37	3.5mm Cannulated Conical Screw 70mm Lng	2
00-2359-012-35	3.5mm Locking Screw 12mm Lng	2
00-2359-014-35	3.5mm Locking Screw 14mm Lng	2
00-2359-016-35	3.5mm Locking Screw 16mm Lng	2
00-2359-018-35	3.5mm Locking Screw 18mm Lng	2
00-2359-020-35	3.5mm Locking Screw 20mm Lng	2
00-2359-022-35	3.5mm Locking Screw 22mm Lng	2
00-2359-024-35	3.5mm Locking Screw 24mm Lng	2
00-2359-026-35	3.5mm Locking Screw 26mm Lng	2
00-2359-028-35	3.5mm Locking Screw 28mm Lng	2
00-2359-030-35	3.5mm Locking Screw 30mm Lng	2
00-2359-032-35	3.5mm Locking Screw 32mm Lng	3
00-2359-034-35	3.5mm Locking Screw 34mm Lng	3
00-2359-036-35	3.5mm Locking Screw 36mm Lng	3
00-2359-038-35	3.5mm Locking Screw 38mm Lng	3
00-2359-040-35	3.5mm Locking Screw 40mm Lng	3
00-2359-042-35	3.5mm Locking Screw 42mm Lng	3
00-2359-044-35	3.5mm Locking Screw 44mm Lng	3
00-2359-046-35	3.5mm Locking Screw 46mm Lng	3
00-2359-048-35	3.5mm Locking Screw 48mm Lng	3
00-2359-050-35	3.5mm Locking Screw 50mm Lng	3
00-2359-052-35	3.5mm Locking Screw 52mm Lng	2
00-2359-054-35	3.5mm Locking Screw 54mm Lng	2
00-2359-056-35	3.5mm Locking Screw 56mm Lng	2
00-2359-058-35	3.5mm Locking Screw 58mm Lng	2
00-2359-060-35	3.5mm Locking Screw 60mm Lng	2
00-2359-065-35	3.5mm Locking Screw 65mm Lng	2
00-2359-070-35	3.5mm Locking Screw 70mm Lng	2
00-2359-010-27	2.7mm Locking Screw 10mm Lng	3
00-2359-012-27	2.7mm Locking Screw 12mm Lng	3
00-2359-014-27	2.7mm Locking Screw 14mm Lng	3
00-2359-016-27	2.7mm Locking Screw 16mm Lng	3
00-2359-018-27	2.7mm Locking Screw 18mm Lng	3
00-2359-020-27	2.7mm Locking Screw 20mm Lng	3
00-2359-022-27	2.7mm Locking Screw 22mm Lng	4
00-2359-024-27	2.7mm Locking Screw 24mm Lng	4
00-2359-026-27	2.7mm Locking Screw 26mm Lng	4
00-2359-028-27	2.7mm Locking Screw 28mm Lng	4
00-2359-030-27	2.7mm Locking Screw 30mm Lng	4
00-2359-032-27	2.7mm Locking Screw 32mm Lng	4
00-2359-034-27	2.7mm Locking Screw 34mm Lng	4
00-2359-036-27	2.7mm Locking Screw 36mm Lng	4
00-2359-038-27	2.7mm Locking Screw 38mm Lng	4
00-2359-040-27	2.7mm Locking Screw 40mm Lng	4
00-2359-042-27	2.7mm Locking Screw 42mm Lng	4
00-2359-044-27	2.7mm Locking Screw 44mm Lng	4
00-2359-046-27	2.7mm Locking Screw 46mm Lng	4
00-2359-048-27	2.7mm Locking Screw 48mm Lng	4
00-2359-050-27	2.7mm Locking Screw 50mm Lng	4

00-2359-052-27	2.7mm Locking Screw 52mm Lng	3
00-2359-054-27	2.7mm Locking Screw 54mm Lng	3
00-2359-056-27	2.7mm Locking Screw 56mm Lng	3
00-2359-058-27	2.7mm Locking Screw 58mm Lng	3
00-2359-060-27	2.7mm Locking Screw 60mm Lng	3

Prod. No.	Description	Set Qty.
00-2359-000-01	5.5mm/4.5mm Locking Screw Set	
00-2359-050-55	5.5mm Cannulated Locking Screw 50mm Lng	4
00-2359-055-55	5.5mm Cannulated Locking Screw 55mm Lng	4
00-2359-060-55	5.5mm Cannulated Locking Screw 60mm Lng	4
00-2359-065-55	5.5mm Cannulated Locking Screw 65mm Lng	4
00-2359-070-55	5.5mm Cannulated Locking Screw 70mm Lng	4
00-2359-075-55	5.5mm Cannulated Locking Screw 75mm Lng	4
00-2359-080-55	5.5mm Cannulated Locking Screw 80mm Lng	4
00-2359-085-55	5.5mm Cannulated Locking Screw 85mm Lng	4
00-2359-090-55	5.5mm Cannulated Locking Screw 90mm Lng	2
00-2359-095-55	5.5mm Cannulated Locking Screw 95mm Lng	2
00-2359-100-55	5.5mm Cannulated Locking Screw 100mm Lng	2
00-2359-050-56	5.5mm Cannulated Conical Screw 50mm Lng	2
00-2359-055-56	5.5mm Cannulated Conical Screw 55mm Lng	2
00-2359-060-56	5.5mm Cannulated Conical Screw 60mm Lng	2
00-2359-065-56	5.5mm Cannulated Conical Screw 65mm Lng	2
00-2359-070-56	5.5mm Cannulated Conical Screw 70mm Lng	2
00-2359-075-56	5.5mm Cannulated Conical Screw 75mm Lng	2
00-2359-080-56	5.5mm Cannulated Conical Screw 80mm Lng	2
00-2359-085-56	5.5mm Cannulated Conical Screw 85mm Lng	2
00-2359-090-56	5.5mm Cannulated Conical Screw 90mm Lng	2
00-2359-012-45	4.5mm Locking Screw 12mm long	2
00-2359-014-45	4.5mm Locking Screw 14mm long	2
00-2359-016-45	4.5mm Locking Screw 16mm long	2
00-2359-018-45	4.5mm Locking Screw 18mm long	2
00-2359-020-45	4.5mm Locking Screw 20mm long	2
00-2359-022-45	4.5mm Locking Screw 22mm long	2
00-2359-024-45	4.5mm Locking Screw 24mm long	2
00-2359-026-45	4.5mm Locking Screw 26mm long	2
00-2359-028-45	4.5mm Locking Screw 28mm long	4
00-2359-030-45	4.5mm Locking Screw 30mm long	4
00-2359-032-45	4.5mm Locking Screw 32mm long	4
00-2359-034-45	4.5mm Locking Screw 34mm long	4
00-2359-036-45	4.5mm Locking Screw 36mm long	4
00-2359-038-45	4.5mm Locking Screw 38mm long	4
00-2359-040-45	4.5mm Locking Screw 40mm long	4
00-2359-042-45	4.5mm Locking Screw 42mm long	4
00-2359-044-45	4.5mm Locking Screw 44mm long	4
00-2359-046-45	4.5mm Locking Screw 46mm long	4
00-2359-048-45	4.5mm Locking Screw 48mm long	4
00-2359-050-45	4.5mm Locking Screw 50mm long	4
00-2359-055-45	4.5mm Locking Screw 55mm long	4
00-2359-060-45	4.5mm Locking Screw 60mm long	2
00-2359-065-45	4.5mm Locking Screw 65mm long	2
00-2359-070-45	4.5mm Locking Screw 70mm long	2



The CE mark is valid only if it is also printed on the product label.



97-2347-109-00 Rev. 1 1004-T08 2ML Printed in USA ©2005, 2007, 2010 Zimmer, Inc.

Contact your Zimmer representative or visit us at www.zimmer.com



zimmer | trauma