



## Zimmer® Periarticular Distal Lateral Femoral Locking Plate



### Plate Features:

- Anatomically contoured plates create a fit that requires little or no additional bending and are helpful for metaphyseal/diaphyseal reduction
- Threaded holes create a 95 degree fixed angle between the plate head and the locking screws to allow screw placement that is parallel to the joint line
- The central locking screw hole in the plate head provides initial reduction of the plate to the condyles
- The low profile plate facilitates fixation without impinging on soft tissues
- Plates are available in a variety of sizes and lengths, from 6 (159mm) to 18 (349mm) holes, left and right
- Dual-compression slots accommodate periarticular screws or conventional stainless steel screws and allow bi-directional compression
- The last diaphyseal plate hole is designed to accommodate the tension device (00-4817-000-05)
- The tapered plate shaft design allows for a minimally invasive technique with submuscular passage of the plate
- Multiple locking holes in the plate head allow placement of the screws to capture fragments
- Thick-to-thin profiles make the plates autocontourable
- The anatomical shape of the head of the plate matches the shape of the distal femur

Prod. No. (Left)	Prod. No. (Right)	Description	Holes	Length
<b>5.5mm Distal Lateral Femoral Locking Plates</b>				
00-2357-102-06	00-2357-101-06	Distal Lateral Femoral Locking Plate	6	159mm
00-2357-102-10	00-2357-101-10	Distal Lateral Femoral Locking Plate	10	223mm
00-2357-102-14	00-2357-101-14	Distal Lateral Femoral Locking Plate	14	287mm
00-2357-102-18	00-2357-101-18	Distal Lateral Femoral Locking Plate	18	349mm

**Designed for use with:**

Zimmer 5.5mm/4.5mm Periarticular Locking Screw Set (00-2359-000-01)  
 Zimmer 5.5mm/4.5mm Periarticular Locking Instrument Set (00-2360-000-01)  
 Zimmer Periarticular Distal Lateral Femoral Plate Standard Jig Set (00-2360-000-11)

**Compatible with:**

Zimmer Periarticular 4.5mm Self-tapping Cortical Screws  
 Zimmer Periarticular 6.5mm Cancellous Screws  
 Zimmer Basic Instrument Set (00-4800-045-20)  
 Zimmer Standard Screw Set (00-4800-045-01)



Prod. No.	Description	Set Qty.
<b>00-2357-000-13 Periarticular Distal Lateral Femoral Locking Plate Set</b>		
00-2357-102-06	Distal Lateral Femoral Locking Plate, 6 Hole, Left	1
00-2357-102-10	Distal Lateral Femoral Locking Plate, 10 Hole, Left	1
00-2357-102-14	Distal Lateral Femoral Locking Plate, 14 Hole, Left	1
00-2357-102-18	Distal Lateral Femoral Locking Plate, 18 Hole, Left	1
00-2357-101-06	Distal Lateral Femoral Locking Plate, 6 Hole, Right	1
00-2357-101-10	Distal Lateral Femoral Locking Plate, 10 Hole, Right	1
00-2357-101-14	Distal Lateral Femoral Locking Plate, 14 Hole, Right	1
00-2357-101-18	Distal Lateral Femoral Locking Plate, 18 Hole, Right	1

Prod. No.	Description	Set Qty.
<b>00-2360-000-11 Periarticular Distal Lateral Femoral Plate Standard Jig Set</b>		
00-2360-090-01	Distal Lateral Femoral Plate Jig, Right	1
00-2360-090-02	Distal Lateral Femoral Plate Jig, Left	1
00-2358-010-00	Distal Lateral Femoral Plate/Jig Case	1

Prod. No.	Description	Set Qty.
<b>00-2360-000-01 5.5mm/4.5mm Periarticular Locking Instrument Set</b>		
00-1147-073-00	Cleaning Stylet	1
00-1147-078-00	Cleaning Brush	1
00-2358-035-00	5.5mm/4.5mm Locking Screw and Instrument Case	1
00-2360-012-01	5.5mm/4.5mm Plate Reduction Instrument	2
00-2360-012-02	5.5mm/4.5mm Plate Reduction Sleeve	2
00-2360-012-03	Plate Reduction Spin Knob	2
00-2360-020-37	3.7mm Standard Cannula	4
00-2360-021-32	3.2mm Standard Cannula	4
00-2360-033-32	3.2mm Standard Drill Tip Guide Wire (5 per box)	2 bx
00-2360-040-45	4.5mm Locking Screw Standard Depth Gauge	1
00-2360-041-55	5.5mm Cannulated Locking Screw Depth Gauge	1
00-2360-053-45	4.5mm Locking Screw Tap	2
00-2360-054-55	5.5mm Cannulated Locking Screw Tap	2
00-2360-065-05	5.0mm Screwdriver Stop Ring	1 bx
00-2360-065-50	5.0mm Hex Standard Screwdriver	2
00-2360-066-50	5.0mm Hex Standard Cannulated Screwdriver	2
00-2360-071-47	4.7mm Standard Cannulated Drill	2
00-2360-080-00	Torque Limiting Attachment	1
00-2360-085-00	Guide Wire Insertor	1
00-2360-086-00	Modular Handle	1
00-2360-087-00	Ratcheting Modular Handle	1
00-2360-088-00	5.5mm/4.5mm Cannula Insertor	1
00-2360-090-04	5.5mm/4.5mm Standard Jig Sleeve	4
00-2360-225-37	3.7mm Standard Drill	2
Also Available:		
00-2360-050-50	5.0mm Hex Forward Captive Screwdriver	0
00-2360-051-50	5.0mm Hex Cannulated Forward Captive Screwdriver	0

Prod. No.	Description	Set Qty.
<b>00-2359-000-01 5.5mm/4.5mm Locking Screw Set</b>		
00-2359-050-55	5.5mm Cannulated Locking Screw 50mm Lng	4
00-2359-055-55	5.5mm Cannulated Locking Screw 55mm Lng	4
00-2359-060-55	5.5mm Cannulated Locking Screw 60mm Lng	4
00-2359-065-55	5.5mm Cannulated Locking Screw 65mm Lng	4
00-2359-070-55	5.5mm Cannulated Locking Screw 70mm Lng	4
00-2359-075-55	5.5mm Cannulated Locking Screw 75mm Lng	4
00-2359-080-55	5.5mm Cannulated Locking Screw 80mm Lng	4
00-2359-085-55	5.5mm Cannulated Locking Screw 85mm Lng	4
00-2359-090-55	5.5mm Cannulated Locking Screw 90mm Lng	2
00-2359-095-55	5.5mm Cannulated Locking Screw 95mm Lng	2
00-2359-100-55	5.5mm Cannulated Locking Screw 100mm Lng	2
00-2359-050-56	5.5mm Cannulated Conical Screw 50mm Lng	2
00-2359-055-56	5.5mm Cannulated Conical Screw 55mm Lng	2
00-2359-060-56	5.5mm Cannulated Conical Screw 60mm Lng	2
00-2359-065-56	5.5mm Cannulated Conical Screw 65mm Lng	2
00-2359-070-56	5.5mm Cannulated Conical Screw 70mm Lng	2
00-2359-075-56	5.5mm Cannulated Conical Screw 75mm Lng	2
00-2359-080-56	5.5mm Cannulated Conical Screw 80mm Lng	2
00-2359-085-56	5.5mm Cannulated Conical Screw 85mm Lng	2
00-2359-090-56	5.5mm Cannulated Conical Screw 90mm Lng	2
00-2359-012-45	4.5mm Locking Screw 12mm long	2
00-2359-014-45	4.5mm Locking Screw 14mm long	2
00-2359-016-45	4.5mm Locking Screw 16mm long	2
00-2359-018-45	4.5mm Locking Screw 18mm long	2
00-2359-020-45	4.5mm Locking Screw 20mm long	2
00-2359-022-45	4.5mm Locking Screw 22mm long	2
00-2359-024-45	4.5mm Locking Screw 24mm long	2
00-2359-026-45	4.5mm Locking Screw 26mm long	2
00-2359-028-45	4.5mm Locking Screw 28mm long	4
00-2359-030-45	4.5mm Locking Screw 30mm long	4
00-2359-032-45	4.5mm Locking Screw 32mm long	4
00-2359-034-45	4.5mm Locking Screw 34mm long	4
00-2359-036-45	4.5mm Locking Screw 36mm long	4
00-2359-038-45	4.5mm Locking Screw 38mm long	4
00-2359-040-45	4.5mm Locking Screw 40mm long	4
00-2359-042-45	4.5mm Locking Screw 42mm long	4
00-2359-044-45	4.5mm Locking Screw 44mm long	4
00-2359-046-45	4.5mm Locking Screw 46mm long	4
00-2359-048-45	4.5mm Locking Screw 48mm long	4
00-2359-050-45	4.5mm Locking Screw 50mm long	4
00-2359-055-45	4.5mm Locking Screw 55mm long	4
00-2359-060-45	4.5mm Locking Screw 60mm long	2
00-2359-065-45	4.5mm Locking Screw 65mm long	2
00-2359-070-45	4.5mm Locking Screw 70mm long	2



0123

97-2347-105-00 10MM Printed in USA ©2005 Zimmer, Inc.

Contact your Zimmer representative or visit us at [www.zimmer.com](http://www.zimmer.com)

